

Life-O-Line
TECHNOLOGIST

FOR COMFORT IN BREATHING

(AN ISO, CE, WHO & GMP CERTIFIED COMPANY)



About us :

The story of Life-O-Line Technologist began over 12 years ago with a few creative individuals working to provide innovative solutions to health care problems. The company was founded in 2008 by Uday Kumar Sinha. Hard work, honesty, knowledge, human face towards patients, opportunity and a bit of luck combined, have led to the creation of a solid company that has consistently provided the products needed for cost-effective, quality health care.

With an innovation history of over 12 years, Life-O-Line Technologist has designed and manufactured products for oxygen delivery and anaesthesia care, Infusion therapy and other equipment for Hospital & ICU departments.

Our efforts are built upon the legacy of product innovation, and to expand the business and create a growing distribution network. In 2012, the company consolidated its manufacturing operations from a few products manufacturing base and mold shops to one state of the art manufacturing plant at a centralized location in Ahmedabad, India. The company has now grown to a dedicated workforce of over 150 employees and it is continuing along its growth and expansion trajectory.

In 2015 the company set up Life-O-line Africa, a distributorship of its quality and dependable products in Africa that is based in South Africa.

CATEGORY INDEX

INFUSION THERAPY

- I V Set Std (Latex Free)
- I V Set Premium (Dehp Free Tubing)
- Blood Transtution Set (Dehp Free Tubing)
- I V Canulla
- Measure Volume Set
- Dial Flow Regulator
- Extension Line
- Pressure Monitoring Line
- Codan Set
- A V Fistulla Needle

ANESTHESIA & RESPIRATORY CARE

- Nasal Cannula
- Oxygen Mask
- Hi-concentration Mask / Adult / Paed
- Nebuliser With T & Mount Pc
- Nebulizer Mask - Adult / Paed
- Nebuliser Kit - Adult
- Multi Flow Ventury Mask Adult
- Single Dial Flow Ventury Mask Adult / Paed
- T - Oxygenator Kit With Tubing
- 3 Ball Siprometr
- Air Cushion Mask
- Bi-pap Mask
- Catheter Mount Std / Exp.
- Hme Filter Adult / paed
- Bvf Filter Adult / Paed
- Hme/bvf Fliter Adult / Paed
- Ventilator Circuit Adult / Paed Neo (Plain)
- Ventilator Circuit Adult / Paed (SWT)
- Ventilator Circuit Adult / Adult (DWT)
- Ventilator Kit Adult / Paed / Neo
- Bain Circuit Adult / Paed
- Guduel Airways (Sizes 000,00,0 To 5)
- Ambu Baug Adult / Paed / Neo
- Endotracheal Tube Plain (sizes 2 To 5.5)
- Endotracheal Tube Cuffed (sizes 6 To 9)
- Reinforced Endotracheal Tube
- Central Venous Catheter (CVC)
- Laryngeal Mask

GASTROENTROLOGY

- Infant Feeding Tube
- Mucus Extractor
- Ryels Tube

UROLOGY & NEPHROLOGY

- TUR Set
- Nelaton Catheter
- Folley's Catheter
- D J Stent
- Urethral Catheter
- Hemodialysis Catheter Kit / Singel /
- Duble/ Triple
- Femoral Catheter
- Blood Tubing Set
- Urometer With Urine Bag
- Urine Bag
- Pcn Catheter
- Pigtail Catheter
- Guide Wire
- Trasnducer Protector
- Nephrostomy Kit

SURGERY & DRAINAGE SYSTEM

- Thoracic Drainage / Trocar Catheter
- Closed Wound Suction Set / Mini Wound
- Suction Set
- Silicon Drainage System
- Yankur Suction Set
- Abdominal Drainage Kit
- Fogarty Catheter
- Water Seal Drainage System / Icd Bag
- Introducer Needle
- Suction Catheter
- Cord Clamp

MISCELLENIOUS


- Surgical Gloves
- Examination Gloves
- Nebuliser Compressor Machine
- Oxygen Flow Meter
- Blood Pressure Instrument
- Biopsy Gun
- Skin Stapler
- Cautery Pencil (Electro Surgical Pencil)
- Epidural Kit
- Spinal Needle
- Bougie
- ID Belt

INFUSION THERAPY




Infusion Set (Standard)

Without Air Vent Spike, Approx 20ml Drops/ml




Infusion Set (Premium)

With Air Vent Spike & Bacteria Retention air inlet with snap on cap "Y" injection port designed for multiple injection, Latex Free.



Blood Transfusion Set

Single Drip Chamber Latex Bulb With "Y" Connector



I V Cannula

(SIZES 14G TO 26G)
With Wings & Injection Port
(18g, 20g, 22g, 24g)




Measure Volume Set (110ml & 150ml)

110ml & 150ml Cylinder With Calibrated Measured Chamber With Blood Marking Scale. Approx 60 Drops/ml Soft & Kink Resistant Tubing.




Dial Flow Regulator

With Or Without I V Set. Flow Control Of I V Fluids With Range Of 5 To 250ml/hrs With "Y" Injection Port With Soft & Kink Resistant Tubing



Extension Line

(SIZES 10CM TO 200CM)
With or Without 3 Way Stop Cock.
With 10cm To 200cm Length Sizes.
For Pressures of 65 PSI.




Pressure Monitoring Line

(SIZES 50CM TO 200CM)
With Male To Male & Male To Female Connector. Kink Resistant Tubing With 50cm To 200cm With 800 PSI.



Codon Set


With 0.2 Micron IV inline filter in order to prevent any particles, entrapped air and bacteria & fungi passing into the venous system, DEHP Free, Latex Free



AV Fistula Needle

Needle To Minimize Interruption For Blood Flow. Stainless Steel Needle With DEHP Free Tubing
(15g, 16g, 17g)

ANESTHESIA & RESPIRATORY CARE



Nasal Cannula (Adult/Pead/Neo)

Kink Resistant Tubing With Star Lumen of 200 cm Length




Oxygen Mask (Adult/Pead)

Kink Resistant Tubing With Star Lumen With Swivel Connector for patient comfort (Adult/Peadiatric)



High Concentration Mask

Kink Resistant Tubing With Star Lumen With Reserver Bag (Adult/Peadiatric)



Nebuliser Mask With T Pcs

With Chamber, Mouth & T Pcs and With Kink Resistan Tubing and also suitable with ventilator circuit for nebulization.



Nebuliser Mask (Adult/Pead)

With Mask, Chamber & Tubing (Adult/Peadiatric)



Nebuliser Kit Adult/Pead)

With Complete Set of Mask, Chamber, Tubing, Mouth Pc , T-pcs With Star Lumen Tubing (Adult/Peadiatric)



Multi Flow Ventury Mask (Adult/Pead)

With Different Oxygen Connectors In Adult & Peadiatric Masks With Kink Resistant Tubing



Multiflow Single Dial Ventury Mask

Single Dial Ventury Mask with kink Resistance Star Lumen Tube with Rotating 3 conn. all % control. Adult / Child



T- Oxygenator

With or Without Tubing & Ventury Valve. With T Piece Of 22mm & 15cm & Suction Port& 40% , 60% Oxygen Conn.



3 Ball Spirometer

Breathing Exerciser With 3 Stage 600cc, 900cc & 1200cc



Air Cushion Mask (Size 0 TO 5)

With PVC Materials & Colour Code Rings (Size 0 To 5)



Bi-Pap Mask (Small/Medium/Large)

With Harness Belt In All Sizes (Small, Medium & Large)

ANESTHESIA & RESPIRATORY CARE



**Catheter Mount
Adult/Pediatric
(Std & Expandable)**
Compatible With ET Tube




**HME Filter Adult /
Pediatric**
Heat Moisture Exchanger With CO₂
Monitoring Port




**BVF Filter Adult /
Pediatric**
With 0.02 Micron Different Layers
of Filters




**HME/ BVF Filter Adult/
Pediatric**
Heat Moisture With Bacteria Viral
Filter With CO₂ Port




**Rebreathing Circuits
Adult/Pediatric/Neonatal**
Corrugated Tubing Of Length Of
1.6 Mtr With Or Without "I"
Connector.



**Rebreathing Circuits
Adult/Pediatric (SWT)**
Corrugated Tubing Of Length Of
1.6 Mtr With Or Without "I"
Connector. With Single Water Trap
Chamber.



**Rebreathing Circuits
Adult/ Pediatric (DWT)**
Corrugated Tubing Of Length Of
1.6 Mtr With Or Without "I"
Connector. With Double Water
Trap Chamber.



**Ventilator Kit - Adult/
Pediatric /Neonatal**
Corrugated Tubing Of Length Of
1.6 Mtr With Or Without "I"
Connector (Plain, SWT & DWT)
with Breathing Filter & Swivel Mount
Catheter.




**Bain Circuit Adult /
Pediatric**
With Rebreathing Bag (2.0 Ltr &
500ml) with Oxygen Tubing &
Corrugated Pipe With ET Tube
Connectors. Adult/pediatric



Guedel Airways
(SIZE 000 TO 5)
Airway Path With Color Coded For
Identification



Ambu Bag - Silicon
Resuscitator, star Lumen
Tubing, reservoir Bag, silicon Mask
& Pressure Valve.
(Adult/Pediatric/Neonatal)



E T Tube Plain
(SIZES 2 TO 5.5)
PVC Tubing With High Volume &
Low Pressure. Sizes Plain 2 To 5.5.
Sizes.

ANESTHESIA & RESPIRATORY CARE



E T Tube Cuffed

(SIZES 5 TO 9)

With Soft Cuff Based On High Volume & Low Pressure. Sizes Plain 2 To 5.5. Sizes Cuffed 5 To 9.



Re-Inforced E.T. Tube

Re-Inforced E.T. Tube - Cuffed & Plain available in various sizes 5.0 to 8.5FG



Central Venous Catheter / Central Line (Single/Double/Triple)

Made up of Specially Formulated Polyurethane material providing strength and longevity during insertion and softens at body temperature.



Laryngeal Mask

(SIZES 1 TO 5)

With PVC Material With 15 MM Connector

GASTROENTROLOGY



Infant Feeding Tube

(SIZES 5FG TO 10FG)

For Neonatal & Peadiatric Feeding



Mucus Extractor

For Collecting Mucus Specimen




Ryels Tube

(SIZES 6FG TO 24FG)


For Nutrition Of Intestinal Secretion

UROLOGY & NEPHROLOGY



TUR Set


For Intigration During TURP Porcess



Nelaton Catheter

(SIZES 6 TO 24FG)


For Short Term Catheterization For Male



Folleys Catheter (2 Way / 3 Way, Latex & Silicon)

SIZES 12 TO 24FG


2Way, 3 Way, Latex and Silicon Foley Balloon Catheter



D J Stent (Closed & Open Tip)

Short & Long Term Connection Of Urethra And Bladder


UROLOGY & NEPHROLOGY



Urethral Catheter
(SIZES 10 TO 14FG)
Short Term Catheterization For Female



Hemodialysis Catheter Kit (Single/Double/Triple)
Made With Polyurethane Material For Strength During Process With Jugular / Straight Tip



Femoral Catheter
Soft tip, Small rotating fixation wing, available in variety of length, size and shapes




Blood Tubing Set
Highly flexible and optimum quality blood tubing set for Nephrology



Urometer With Urine bag
Urine Bag with Urometer designed for urine measurement on hourly basis during examination, available with sampling port and "T" outlet Valve.




Urine Bag (Super/Premium/Deluxe)
With Non Return Valve & 2000ml Capacity



PCN Catheter
SIZES 08, 10 & 12F
Percutaneous Nephrostomy Catheter available with our without needle. Durable with superior bio-compatibility.



Pigtail Catheter
SIZES 08, 10 & 12F
Highly Radiopaque For Better Visualization Tapered at Both Ends For Easy Insertion



Guide Wire (PTFE/Nitinol)
PTFE Coated or Nitinol Guide wire available in various sizes



Transducer Protector
Durable and Fine-Finish Transducer protector for Haemodialysis.



Nephrostomy Kit
Nephrostomy Set Used For Percutaneous Placement Of A Pigtail Or Malecot Catheter In The Renal Pelvis For Nephrostomy Drainage

SURGERY & DRAINAGE SYSTEM



Thoracic Drainage / Trocar Catheter

(SIZES 8 TO 40)

Trocar With Catheter For Placement & Handling



Closed Wound Suction Set / Mini Wound Suction Set

(SIZES 6 TO 18)


Universal Trocar And X-ray Opaque Line On The Catheter. Also Available In Mini Size For Paediatric Use.



Silicone Drainage System

(SIZES 8 TO 40)

For Drainage After Thoracic Surgery



Yankauer Suction Set (Plain Tip/Crown Tip)


Handle With Siliconised Tubing For Suction Of Fluids



Abdominal Drainage Kit

(SIZES 20 TO 36)

Soft Catheter With Collection Bag Capacity Of 2000ml



Fogarty Catheter

Embolectomy Catheter, Fogarty with 4 Fr. x 80 Cm, 9 mm Balloon Size Arterial PVC Body/Latex Balloon



Water Seal Drainage System / ICD Bag


SIZE 1000ML & 2000ML

Specially Designed Drainage Bag Suitable For Pleural Drainage In Conjunction With Intercoastal Drainage Catheters



Introducer Needle

Extremely Low Insertion Forces Increase Accessibility And Reduce Tissue Trauma. Tapered Hub Inner Diameter Facilitates Insertion Of The Guide Wire Smoothly



Suction Catheter

SIZE 6.0 TO 20.0FG

Designed For Tracheal / Bronchial Suction. Specially Designed For Easy Passage Through Endotracheal and Tracheotomy Tubes With Minimum Irritation To Patient



Cord Clamp

Clamping For Umbilical Cord

MISCELLENOUS



Surgical Gloves

(SIZES 6 TO 8)

Made up of high quality natural latex. Powdered or powder-free available as per requirement



Examination Gloves (Small/Medium/Large)

Medium weight examination gloves for medical and lab purpose



Nebuliser Compressor Machine

Compressor With Disposable Side Stream, Mouthpiece, Tubing, Filter And Mask (adult / Paed.)




Oxygen Flow Meter

Back pressure compensated oxygen flow meter with accurate gas flow measurement.



Blood Pressure Instrument

300 mm / Hg Blood Pressure Equipment For Clinical Use



Biopsy Gun

Single Use, Guillotine-type, Soft Tissue Biopsy Instrument Which Performs The Entire Bioptic Manoeuvre With Precision.



Skin Stapler

Special Designed and Manufactured Disposable Skin Stapler For Surgical Skin Suture



Cautery Pencil

3 Mtrs. Length Electrosurgical Pencil



Epidural Kit

(SIZES 16 & 18NG)


Epidural Catheter Kit Designed For Short Term & Long Term Continuous Anaesthesia & Pain Therapy



Spinal Needle


(SIZES 18 TO 27G)

Easy To Dispose Off Needles



Bougie

Disposable Bougie Useful in Processes Involving Post Extubation as Well As Reintroduction of Tubes



ID Belt

Designed For Managing Moderate To Heavy Incontinence In Active And Independent People

MD-5 (A & B Category)

FORM MD-5

[Please submit (6) of rule 20 and sub-rule (6) of rule 28]

Licence to Manufacture for Sale or for Distribution of Class A or Class B medical device

Licence Number: MF/CD/01/2018/000001

1. Mfr Life O-Line Technology, Plot No. 864/1, Pank, Nr. ICCL Petrol Pump, Ahmedabad-Mehmedabad Road/PJ And AT Hirapur Village, Ahmedabad, Gujarat (India) - 382445 Telephone No. 9898162376 FAX: 9898162376

2. Details of medical device(s) (Annexure)

3. This licence is subject to the provisions of the Medical Device Rules, 2017 and conditions prescribed therein.

S.No. Details Of Device(s)

1. Generators/Pressure MONITORING LINE Model No. PM-25PM-50PM-100PM-250 - NA Intended Use: Monitoring a fluid intravenous infusion Treatment Class of medical device Class B Material of construction PVC, PC, HDPE Dimension(s) any/NA Shelf life 3 Years Details or Non sterile/ Sterile Brand Name(s) registered under the Trade Marks Act, 1999: Life-O-Line Generators/Pressure MONITORING LINE WITH 1 WAY STOP COCK Model No. ECT/100ECT/250ECT/50ECT/100ECT/250ECT - NA Intended Use: The intended use of the Short Term Catheter Catheter allowing temporary/Short Term (less than 30 days) vascular access for dialysis Class of medical device Class C Material of construction PVD, PPS/HDPE Dimension(s) any/NA Shelf life 3 Years Details or Non sterile/ Sterile Brand Name(s) registered under the Trade Marks Act, 1999: Life-O-Line

MD-9 (C & D Category)

FORM MD-9

[Please submit (7) of rule 28]

Licence to Manufacture for Sale or for Distribution of Class C or Class D medical device

Licence Number: MF/CD/01/2018/000012

1. Mfr Life O-Line Technology, Plot No. 864/1, Pank, Nr. ICCL Petrol Pump, Ahmedabad-Mehmedabad Road/PJ And AT Hirapur Village, Ahmedabad, Gujarat (India) - 382445 Telephone No. 9898162376 FAX: 9898162376

2. Details of medical device(s) (Annexure)

3. The names, qualifications and experience of the competent technical staff responsible for the manufacture and testing of the above mentioned medical device(s) As per records maintain by the manufacturer

4. This licence is subject to the provisions of the Medical Device Rules, 2017 and conditions prescribed therein.

S.No. Details Of Device(s)

1. Generators/DOUBLE LUMEN DIALYSIS KIT / JUGULAR / STRAIGHT Model No. DL-1213DL-1215DL-1215DL-1220 - NA Intended Use: The intended use of the Short Term Catheter Catheter allowing temporary/Short Term (less than 30 days) vascular access for dialysis Class of medical device Class C Material of construction PVD, PPS/HDPE Dimension(s) any/NA Shelf life 3 Years Details or Non sterile/ Sterile Brand Name(s) registered under the Trade Marks Act, 1999: Life-O-Line Generators/DOUBLE LUMEN DIALYSIS KIT Model No. JUG-1212JUG-1215 - NA Intended Use: The intended use of the Short Term Catheter Catheter allowing temporary/Short Term (less than 30 days) vascular access for dialysis Class of medical device Class C Material of construction PVC/PPS/HDPE Dimension(s) any/NA Shelf life 3 Years Details or Non sterile/ Sterile Brand Name(s) registered under the Trade Marks Act, 1999: Life-O-Line

ISO 9001 : 2015

Certificate of Registration

This is to certify that the Quality Management System of

Life-O-Line Technology

at

864/1, Nr. Indian Petrol Pump, Hirapur Char Rasta, Mehmedabad Road, Ahmedabad Gujarat, 382435, India

has been independently assessed and is compliant with the requirements of:

ISO 9001:2015

For the following scope of activities:

Manufacturing and Marketing of Medical Disposables, Respiratory Care Products and Equipments

Certificate Number: UQ-12311

Validity of this certificate can be verified at www.ukcertifications.org.uk/verify

Date of Initial Registration 27th December 2018

Date of this Certificate 27th December 2018

Certificate Expiry 26th December 2021

Recertification Due subject to the company maintaining its system to the required standard 26th December 2021

Authorised Signatory

ISO 13485 : 2016

Certificate of Registration

This is to certify that the Quality Management System for Medical Device of

Life-O-Line Technology

at

864/1, Nr. Indian Petrol Pump, Hirapur Char Rasta, Mehmedabad Road, Ahmedabad Gujarat, 382435, India

has been independently assessed and is compliant with the requirements of:

ISO 13485:2016

For the following scope of activities:

Manufacturing and Marketing of Medical Disposables, Respiratory Care Products and Equipments

Certificate Number: UQ-12310

Validity of this certificate can be verified at www.ukcertifications.org.uk/verify

Date of Initial Registration 27th December 2018

Date of this Certificate 27th December 2018

Certificate Expiry 26th December 2021

Recertification Due subject to the company maintaining its system to the required standard 26th December 2021

Authorised Signatory

CE

Certificate of Compliance

CE

Confirms that the product: Medical Disposable & Respiratory Care Products of class I according to the Medical Council Directive 93/42/EEC of June 1993

Certificate No.: CE-3739

Manufacturer : Life-O-Line Technology

Address : 864/1, Nr. Indian Petrol Pump, Hirapur Char Rasta, Mehmedabad Road, Ahmedabad Gujarat, 382435, India

Product Description : As Per Annexure 1

Complies with the requirements applicable to it

The manufacturer's technical documentation as required for class I device has been reviewed and found to comply with the requirements any significant changes in the design or construction of the product, not agreed upon by us. This declaration will lose its validity.

This certificate is based upon the following conditions:

1. It applies only to the quality system maintained in the manufacture of above mentioned product and is not intended for design or implementation purposes of regulatory compliance.

2. The certificate remains valid until the manufacturing conditions or the quality system are changed.

3. The certificate validity is confirmed by positive results or surveillance audits.

4. After failing the required CE regulation, the manufacturer shall inform us within 10 days of the date of failure.

5. The CE mark on design above can be used under the responsibility of the manufacturer after declaration of an EC Declaration of conformity and compliance with all relevant CE Directives. The statement is based on a single evaluation of one sample of above mentioned product. It does not imply an endorsement of the design or construction.

Validity of this certificate can be verified at www.ukcertifications.org.uk/verify

Date of initial registration 27th December 2018

Date of this certificate 27th December 2018

Certificate expiry 26th December 2021

Recertification due subject to the company maintaining its system to the required standard 26th December 2021

Authorised Signatory

WHO / GMP

Certificate of Compliance

Certificate Number: UQ-5531

This is to certify that

Life-O-Line Technology

at

864/1, Nr. Indian Petrol Pump, Hirapur Char Rasta, Mehmedabad Road, Ahmedabad Gujarat, 382435, India

Has successfully implemented the Quality Management System and been found working satisfactorily as per the terms of "Good Manufacturing Practice" as laid down by "World Health Organization" which has been in accordance to the requirements of:

WHO-GMP

Scope: Manufacturing & Marketing of Medical Disposable & Respiratory Care Products

This certificate is issued under the following conditions:

1. It applies only to the quality system maintained in the manufacture of above mentioned Medical Products.

2. The certificate remains valid until the manufacturing conditions or the quality system are changed and is subject to continuous surveillance according to the WHO-GMP Guidelines.

3. The certificate validity is confirmed by positive results or surveillance audits.

Validity of this certificate can be verified at www.ukcertifications.org.uk/verify

Date of Certification 02nd June 2017

1st Surveillance Audit Due 01st June 2018

2nd Surveillance Audit Due 01st June 2019

Certificate Expiry subject to the company maintaining its system to the required standard 01st June 2020

Authorised Signatory





Life-O-Line
TECHNOLOGIST

FOR COMFORT IN BREATHING

(AN ISO, CE, WHO & GMP CERTIFIED COMPANY)

MFG., IMPORTER & EXPORTER OF : MEDICAL DISPOSABLES & MEDICAL DEVICES

Plot No. 864/1, Near Indian Petrol Pump, Hirapur Char Rasta,
Mahemdabad Road, Ahmedabad - 382 435.

E-mail : lifeoline2011@yahoo.com Web : www.lifeoline.com

Mobile : +91 9898 162576 Customer Care : +91 76000 20902